

Messlöffel mit Graduierung - Polypropylen

Artikel-Nummer	281893 281894 281895 281896
Material	Polypropylen
Farbe	naturfarbig
Inhalt	281893 – 5 ml 281894 – 10 ml 281895 – 25 ml 281896 – 50 ml
Gesamtlänge	281893 – 65 mm 281894 – 70 mm 281895 – 85 mm 281896 – 100 mm
Teilung ml	281893 – 1/2,5/5 281894 – 2,5/5/6,5/7,5/8,5/10 281895 – 5/10/15/20/25 281896 – 10/20/30/40/50

