

KERN OZL467



The flexible, affordable all-rounder with zoom function for schools, training companies, inspection authorities and laboratories



Piktogramme

STANDARD



3 YEARS
WARRANTY

Power supply

Input voltage:	100 - 240 V
Mains adapter/adapter included:	EURO

Category

Brand:	Optics
--------	--------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	460 x 410 x 290 mm
Gross weight:	5,500 kg
Net weight:	4,585 kg
Shipping method:	Parcel service

Construction

Dimensions (WxDxH):	300 x 240 x 420 mm
Eye distance (max):	75
Eye distance (min):	55
Optical system:	Greenough
Rotation of tube 360°:	yes
Tube:	Binocular
Tube angle:	45°