

**PRODUCT SPECIFICATION**

**Bowls grey, Polypropylene, 0,5 to 11,0 liters**

Item numer	309650 – 0,5 liters, 309651 – 1,0 liters 309652 – 1,5 liters, 309653 – 2,5 liters 309654 – 5,0 liters, 309655 – 8,0 liters 309656 – 11,0 liters
nominal volume	500 to 11000 ml
material	Polypropylene
colour	Granite gray
Weight/pcs.	309650 – 0,028 kg, 309651 – 0,044 kg 309652 – 0,071 kg, 309653 – 0,109 kg 309654 – 0,207 kg, 309655 – 0,247 kg 309656 – 0,351 kg
Height	309650 – 75 mm, 309651 – 85 mm 309652 – 100 mm, 309653 – 110 mm 309654 – 120 mm, 309655 – 130 mm 309656 – 140 mm
Diameter	309650 – 130 mm, 309651 – 160 mm 309652 – 200 mm, 309653 – 230 mm 309654 – 300 mm, 309655 – 350 mm 309656 – 400 mm
quantity per pack	0,5 to 2,5 liters – 24 pcs 5,0 to 11,0 liters – 12 pcs

